

EQUIPMENT OVERVIEW 2011

Summary: iLs offers both Receptive and Interactive Language systems. All receptive systems – the Pro, Total Focus, Focus 200 and Kick Start - come complete with auditory and movement equipment. The Interactive Language equipment is an add-on product which is primarily used by clinicians helping those with auditory processing challenges.

iLs Pro is a portable air/bone conduction system which requires completion of the 4-day Advanced Certification course and required case studies before purchasing.

- Price: \$1950
- Full library of processed music loaded onto a 160 GB iPod for customizing receptive programs
- Custom, rechargeable amplifier – 7+ hours of continual use
- Sennheiser air/bone conduction headphones
- Integration Kit (movement & visual equipment/guidebook)
- Waist pack for portability; iPod and amp rechargers



The Pro differs from the Total Focus and Focus 200 systems in three ways: 1- the Pro has greater customization options, higher bone conduction level options and left/right ear balance control.

iLs Total Focus – premier air/bone conduction system for use in clinic, hospital, school or home

- ❖ Wholesale price - \$1395 Suggested Retail price - \$1795
- ❖ Pre-loaded Programs:
 - 1) Sensory Motor
 - 2) Concentration/Attention
 - 3) Reading/Auditory Processing
 - 4) Optimal Performance for teens/adults
 - 5) Interactive Language Activities
- ❖ Equipment includes: iPod loaded with 5 complete programs (over 200 hours of processed music), mini-amplifier to power the bone conduction headphones; Sennheiser headphones retro-fitted with bone conductor and high quality cabling; adjustable waist pack for all ages; Playbook and Integration Kit



iLs Focus 200 – same product as the Total Focus with only 2 of the 5 programs

- ❖ Wholesale price \$1195 Suggested Retail price \$1495
- ❖ Pre-loaded Programs: choose any 2 from the 5 Total Focus programs listed above
- ❖ Equipment includes: iPod loaded with 2 complete programs, mini-amplifier to power the bone conduction headphones; Sennheiser headphones retro-fitted with bone conductor and high quality cabling; adjustable waist pack; Playbook and Integration Kit

iLs Kick Start is designed for home use as a follow-up, or supplement, to clinic therapy. It is a portable system which has air conduction only.

- ❖ Wholesale Price: \$500 Retail: \$699
- ❖ 12-week program which can be customized per user needs
- ❖ 4 GB Nano with armband
- ❖ Sennheiser air conduction headphones
- ❖ Integration Kit (equipment/guidebook for visual tracking & movement)



Integration Kit - The Pro, Total Focus, Focus 200 and Kick Start include all equipment necessary for iLs' multi-sensory program: auditory equipment and Integration Kit (pictured right) - balls, bean bags, balance board & Playbook guide for visual/vestibular activities.



The iLs Interactive Language Program is designed to help children and adults with various auditory processing skills such as **phonemic awareness, auditory memory, auditory figure ground and dichotic listening**.

Interactive Language training is usually incorporated into the latter half of iLs programs. It can be implemented by any trained teacher or professional. The Interactive Language activities are played from the Focus or Pro iPods, or from a dedicated Nano used expressly for the activities.



Interactive Language Kit – microphone, stand, amplifier, Nano & docking station - \$625

Interactive Language Kit w/o Nano – microphone, stand, amplifier, docking station w/o Nano - \$425

Interactive Language Exercises

1. **Children's Songs:** getting to know one's own voice
2. **Vocal Production:** singing, humming, toning
3. **Repeating Words:** focus on enunciation
4. **Sibilants:** high frequency sounds in words such as 's' and 'ch' (fricatives)
5. **Repeating Phrases:** focus on auditory memory
6. **Tongue Twisters:** enunciation and fun
7. **Auditory Memory:** 2, 3, 4 and 5-word lists
8. **Filtered Words:** filtering out lower frequencies to improve auditory focus
9. **Auditory Figure Ground:** words/phrases with competing background noise
10. **Dichotic Words:** left and right ear perception and processing